

Radioactivity

- 8.1- Scintillation counter, Beckman Coulter, LS 6500
- 8.2- Geiger-Müller counter, Ludlum Measurements, Inc., Model 3
- 8.3- Table-top centrifuge, Beckman, Avanti J-15R
- 8.4- Hybridization oven, single chamber, Techne, FHB2DP

The Radioactivity Room is located at ABI 359. ABI will provide the annual recertification of the scintillation counter.

Revision #3

Created 2025-09-12 18:23:36 UTC by Rodney Harris

Updated 2026-01-06 20:26:32 UTC by Rodney Harris